

Yellow Latex Powder Free Lemon Scent Powder Free Natural Rubber Latex Examination Glove



Benefits

- High visibility due to yellow color
- Pleasant, fresh Lemon Scent
- Tactile sensation like second skin
- Excellent elasticity
- Outstanding durability & stress relaxation
- Moderate texture

Features

- Low protein content <50µg/dm²
- AQL 1.5
- Product is tested according to ASTM & EN testing methods for examination gloves and meet ASTM D3578 & EN 455 1;2;3 international standard requirements

Powder Free: yes **Surface:** Textured **Cuff:** Beaded

Colors:



| Size | Palm Width (mm) | Length (mm) |
|-------------|-----------------|-------------|
| Extra Small | < 80 | Min. 240 |
| Small | 80 +/- 10 | Min. 240 |
| Medium | 95 +/- 10 | Min. 240 |
| Large | 110 +/- 10 | Min. 240 |
| Extra Large | > 110 | Min. 240 |

| Physical Properties | Before ageing | After ageing | Glove Thickness | (Single Wall/mm) |
|---------------------|---------------|--------------|-----------------|------------------|
| Tensile Strength | Min. 18 MPa | Min. 14 MPa | Cuff | Min. 0.08 |
| Elongation at break | Min. 650% | Min. 500% | Palm | Min. 0.11 |
| Force at break | Min. 6.0 N | Min. 6.0 N | Finger | Min. 0.12 |

| Inspection Parameter | Defect Description | Inspection Plan | Level AQL |
|--------------------------|--|-----------------|--------------|
| Dimension | Glove Thickness (Cuff, Palm, Fingers, Width, Length) | | S-2 4.0 |
| Physical Properties | Tensile Strength, Ultimate Elongation Force at Brake | | S-2 4.0 |
| Visual Inspection | Major Defect, Minor defects | | G-1 2.5, 4.0 |
| 1000 ml Water Tight Test | Holes, Leaks | | G-1 1.5 |
| Powder Free Residue | Powder residue content | | N = 5 |
| Protein | Protein Content | | N = 3 |

Packing Description

| | |
|------------------------|--|
| Packing | 100 pieces per Dispenser, 10 x 100 pieces per Carton |
| Lot Number Structure | (YYYY = Year of manufacturing, XXX = Julian Day) |
| Recommended Shelf Life | 3 Years upon manufacturing date |