

## Purple Nitrile Powder Free 285mm 7MIL Powder Free Nitrile Glove



### Benefits

- Silky comfort fitting
- Soft & strong
- Elastic & yields to hand
- Superior puncture resistance
- Fit with wet & dry grip
- Fingertips textured

### Features

- AQL 1.5
- Product is tested according to ASTM & EN testing methods for examination gloves and meet ASTM D6319 & EN 455 1-4, EN374 (CAT III, Type B) international standard requirements
- Food contact FDA & EN 1186

**Powder Free:** yes    **Surface:** Textured    **Cuff:** Beaded

**Colors:**



Size	Palm Width (mm)	Length (mm)
Small	83 - 86	Min. 285
Medium	93 - 96	Min. 285
Large	102 - 105	Min. 285
XLarge	110-115	Min. 285
XXLarge	118-123	Min. 285

Physical Properties	Before ageing	After ageing	Glove Thickness	(Single Wall/mm)
Tensile Strength	Min. 17 MPa	Min. 16 MPa	Cuff	Min. 0.10
Elongation at break	500-600%	400-500%	Palm	Min. 0.12
Force at break	6.0 N	6.0 N	Finger	Min. 0.17

Inspection Parameter	Defect Description	Inspection Plan	Level	AQL
Dimension	Glove Thickness (Cuff, Palm, Fingers, Width, Length)		S-2	4.0
Physical Properties	Tensile Strength, Ultimate Elongation Force at Brake		S-2	4.0
Visual Inspection	Major Defect, Minor defects		G-1	2.5, 4.0
1000 ml Water Tight Test	Holes, Leaks		G-1	1.5
Powder Free Residue	Powder residue content		N = 5	

### Packing Description

Packing	50 pieces per Dispenser, 10 x 50 pieces per Carton
Lot Number Structure	(YYYY = Year of manufacturing, XXX = Julian Day)
Recommended Shelf Life	3 Years upon manufacturing date