

# Latex Powdered Super Light Weight Powdered Natural Rubber Latex Glove



### Benefits

- Second skin feeling due to thinness & light weight
- Excellent for Food Industry & Retail
- Good elasticity
- Smooth grip

### Features

- Protein content 160 - 200µg/dm<sup>2</sup>
- AQL 4.0
- Product is tested according to international standard requirements

**Powder Free:** yes

**Surface:** Smooth

**Cuff:** Beaded

**Colors:**



Size	Palm Width (mm)	Length (mm)
Extra Small	74-76	Min. 240
Small	82-84	Min. 240
Medium	93-95	Min. 240
Large	103-105	Min. 240
Extra Large	112-114	Min. 240

Physical Properties	Before ageing	After ageing	Glove Thickness	(Single Wall/mm)
Tensile Strength	18 MPa	14 MPa	Cuff	Min. 0.06
Elongation at break	650%	600%	Palm	Min. 0.08
Force at break	< 6 N	< 6 N	Finger	Min. 0.08

Inspection Parameter	Defect Description	Inspection Plan	Level	AQL
Dimension	Glove Thickness (Cuff, Palm, Fingers, Width, Length)		S-2	4.0
Physical Properties	Tensile Strength, Ultimate Elongation Force at Brake		S-2	4.0
Visual Inspection	Major Defect, Minor defects		G-1	4.0, 4.0
1000 ml Water Tight Test	Holes, Leaks		G-1	4.0
Powder Free Residue	Powder residue content		N = 5	
Protein	Protein Content		N = 3	

### Packing Description

Packing	100 pieces per Dispenser, 10 x 100 pieces per Carton
Lot Number Structure	(YYYY = Year of manufacturing, XXX = Julian Day)
Recommended Shelf Life	5 Years upon manufacturing date

Specifications and products shown are correct at time of printing and may be subject to change without prior notice

1210V00\_01