

Nitrile Powder Free Super Light Weight Powder Free Nitrile Examination Glove



Benefits

- Silky comfort fitting
- Soft & light (3.0 g)
- Elastic & yields to hand
- Fit with wet & dry grip
- Fingertips textured

Features

- AQL 1.5
- Chlorinated
- Product is tested according to ASTM testing methods for examination gloves and meet ASTM D6319 & EN 455 1-4, EN374 (CAT III, Type B) international standard requirements

Powder Free: yes **Surface:** Textured **Cuff:** Beaded

Colors:



Size	Palm Width (mm)	Length (mm)
Extra Small	75 ± 5	Min. 240
Small	85 ± 5	Min. 240
Medium	95 ± 5	Min. 240
Large	106 ± 5	Min. 240
Extra Large	116 ± 5	Min. 240
Extra Extra Large	125 ± 5	Min. 240

Physical Properties	Before ageing	After ageing	Glove Thickness	(Single Wall/mm)
Tensile Strength	Min. 18 MPa	Min. 14 MPa	Cuff	Min. 0.04
Elongation at break	Min. 500%	Min. 400%	Palm	Min. 0.05
Force at break	6.0 N	6.0 N	Finger	Min. 0.06

Inspection Parameter	Defect Description	Inspection Plan	Level	AQL
Dimension	Glove Thickness (Cuff, Palm, Fingers, Width, Length)		S-2	4.0
Physical Properties	Tensile Strength, Ultimate Elongation Force at Brake		S-2	4.0
Visual Inspection	Major Defect, Minor defects		G-1	2.5
1000 ml Water Tight Test	Holes, Leaks		G-1	1.5
Residual Powder Level	Residual Powder Level (mg/glove)		N = 5	

Packing Description

Packing	100 or 200 pieces per Dispenser, 10 x 100 or 200 pieces per Carton
Lot Number Structure	(YYYY = Year of manufacturing, XXX = Julian Day)
Recommended Shelf Life	3 Years upon manufacturing date