

## Orange Nitrile 270mm 6mil - High Risk Powder Free Nitrile Examination Glove



### Benefits

- Extra long & Strong
- High visibility
- Chemotherapy Drug Tested
- Silky comfort fitting
- Elastic & yields to hand
- Superior puncture resistance
- Fit with wet & dry grip

### Features

- AQL 1.5
- Product is tested according to ASTM testing methods for examination gloves and meet ASTM D6319 & EN 455 1-4, EN374 (CAT III, Type B), EN1186 international standard requirements

**Powder Free:** yes    **Surface:** Textured    **Cuff:** Beaded

**Colors:**



Size	Palm Width (mm)	Length (mm)
Small	85 +/- 5	Min. 270
Medium	95 +/- 5	Min. 270
Large	106 +/- 5	Min. 270
XLarge	116 +/- 5	Min. 270
XXLarge	125 +/- 5	Min. 270

### Physical Properties

	Before ageing	After ageing
Tensile Strength	Min. 18 MPa	Min 14 MPa
Elongation at break	Min. 500%	Min. 400%
Force at break	Min. 6.0 N	Min. 6.0 N

### Glove Thickness

	(Single Wall/mm)
Cuff	Min. 0.08
Palm	Min. 0.11
Finger	Min. 0.17

### Inspection Parameter

Inspection Parameter	Defect Description	Inspection Plan	Level	AQL
Dimension	Glove Thickness (Cuff, Palm, Fingers, Width, Length)		S-2	4.0
Physical Properties	Tensile Strength, Ultimate Elongation Force at Brake		S-2	4.0
Visual Inspection	Major Defect, Minor defects		G-1	2.5, 4.0
1000 ml Water Tight Test	Holes, Leaks		G-1	1.5
Powder Free Residue	Powder residue content		N = 5	

### Packing Description

Packing	50 pieces per Dispenser, 10 x 50 pieces per Carton
Lot Number Structure	(YYYY = Year of manufacturing, XXX = Julian Day)
Recommended Shelf Life	3 Years upon manufacturing date