

300mm Nitrile - LOW DERMA

Powder Free Nitrile Glove Accelerator Free



Benefits

- Formulated without Chemical Accelerators and Sulphur
- Protect Against Glove Induced Type I & Type IV Hypersensitivities and Blood-borne Pathogens & Viruses
- Maximum comfort & good fit

Features

- AQL 1.5
- Accelerator Free & Low Dermatitis
- Product is tested according to ASTM testing methods for examination gloves and meet ASTM D6319 & EN 455 1:3, EN374 (CAT III, Type B), EN1186 international standard requirements

Powder Free: yes **Surface:** Textured **Cuff:** Beaded

Colors:



Size	Palm Width (mm)	Length (mm)
Extra Small	75 +/- 5	Min. 295
Small	85 +/- 5	Min. 295
Medium	95 +/- 5	Min. 295
Large	106 +/- 5	Min. 295
Extra Large	116 +/- 5	Min. 295

Physical Properties	Before ageing	After ageing	Glove Thickness	(Single Wall/mm)
Tensile Strength	Min. 18 MPa	Min 14 MPa	Cuff	Min. 0.05
Elongation at break	Min. 500%	Min. 400%	Palm	Min. 0.06
Force at break	Min. 6.0 N	Min. 6.0 N	Finger	Min. 0.08

Inspection Parameter	Defect Description	Inspection Plan	Level	AQL
Dimension	Glove Thickness (Cuff, Palm, Fingers, Width, Length)		S-2	4.0
Physical Properties	Tensile Strength, Ultimate Elongation Force at Brake		S-2	4.0
Visual Inspection	Major Defect, Minor defects		G-1	2.5, 4.0
1000 ml Water Tight Test	Holes, Leaks		G-1	1.5
Powder Free Residue	Powder residue content		N = 5	

Packing Description

Packing	100 pieces per Dispenser, 10 x 100 pieces per Carton
Lot Number Structure	(YYYY = Year of manufacturing, XXX = Julian Day)
Recommended Shelf Life	3 Years upon manufacturing date